

3DCultr Lung Cancer Organoid Growth Medium(Human)

Product description

3DCultr Lung Cancer Organoid Growth Medium (Human) is a serum-free medium that can be used for the establishment and long-term culture of lung cancer organoids derived from cells or tissues. In the presence of extracellular matrix, The unique components and abundant cytokines contained in the culture medium can promote the rapid growth and formation of lung cancer organoids from lung cancer cells. The organoid formation process is smooth and rapid, while maintaining high characteristics and vitality of lung cancer cells, which provides a basis for subsequent support for physiological functions, disease research and precision medicine of lung cancer organoids.

Specifications

Catalog Number	C231110E/C231110S/C231110M
Specifications	50 mL/100 mL/500 mL

Components

Component Number	Component Name	C231110E	C231110S	C231110M
C231110-A	Lung Cancer Organoid Growth Medium	45 mL	90 mL	450 mL
С231110-В	Nutritional components (10 $ imes$)	5 mL	10 mL	50 mL

Storage

Stored at -25°C~-15°C, valid for 1 year; when stored at 2~8°C, valid for 1 month.

Notes

1. The operations such as packaging and use of the product should be carried out in a sterile environment, and the experimental equipment (such as: pipette tips, product tubes, etc.) in contact with the product should be pre-cooled before use.

- 2. For your safety and health, please wear a lab coat and disposable gloves.
- 3. For research use only.

Instructions

1. Preparing human lung cancer organoid growth medium

Complete lung cancer organoid culture medium was prepared under sterile operating conditions. The following is the procedure for preparing 100 mL of complete culture medium. If the required 1



amount is different, the amount can be adjusted accordingly.

1) Thaw component B at room temperature or slowly thaw at 2~8°C overnight. Avoid repeated freezing and thawing, prepare and thaw immediately;

2) Take 90 mL of component A out of the refrigerator and return it to room temperature;

3) Add 10 mL of component B to the component A and mix evenly; if not used temporarily, store at 2~8°C for a short period of time.

4) You can add 1% double antibody when using.

2. Primary culture of human intestinal cancer

 Material Collection: After the specimen is removed from the body, collect the material as soon as possible. Use sterile instruments to ensure a sterile environment, place the tumor tissue into a 15 mL centrifuge tube containing 5 mL of primary tissue preservation solution, and transport it at 4°C.

2) Washing: In a biosafety cabinet, remove the sample tube and discard the tissue preservation solution. Add an appropriate amount of cold PBS containing antibiotics, and wash the tissue repeatedly. Discard the PBS.

3) Repeat Washing: Repeat step 2 three times.

4) Tissue Processing: After removing the PBS buffer, move the tissue block to a 10 cm sterile petri dish containing 10 mL of cold primary tissue preservation solution, and cut the tissue into pieces with sterile ophthalmic microscissors (approximately 0.5 mm~1 mm in diameter).

5) Repeat Washing: Wash the tissue fragments with room temperature PBS, repeating step 2 three times.

6) Collect tissue fragments, add tissue digestion solution for digestion for 20~30 minutes, pipe repeatedly and pass through a 70 μ m mesh to collect lung cancer cells. If there are few cells, repeat once.

7) Red Blood Cell Lysis: Add 10 mL of red blood cell lysis buffer and shake on a rocker shaker at room temperature for 10 minutes.

8) Repeat Cleaning: After lysis is completed, use DMEM/F12 at room temperature and repeat step 2 three times.

9) Organoid Seeding Plate: Adjust the cell density to $2 \sim 3 \times 10^6$, mix evenly with Arcegel Matrigel 1:1, seed the cell suspension in a 24-well plate at 40~60 µL per well, and place it at 37°C for 15~30 min, add preheated organoid culture medium, 750 µL to each well.

10) Organoid Culture: Place the culture plate in a 37° C CO₂ incubator. Change the culture medium every 2 days. When adding culture medium, keep the tip facing the side wall and add slowly.

11) Organoid Observation: Observe the organoids and take pictures every day to understand the initial number of organoids, proliferation rate, morphology, microbial contamination, etc.