



TOP10 Chemically Competent Cell

Product description

TOP10 Chemically Competent Cell from YEASEN is the competent cells obtained by the special process of *Escherichia coli* strain TOP10 with Streptomycin (StrR) resistance. It can be used for the chemical transformation of DNA with the conversion efficiency up to 10^8 cfu/ μ g plasmid DNA. It is widely used for efficient DNA cloning and plasmid amplification, and can ensure stable inheritance of high-copy plasmids. At the same time, it can be used to construct clones, blue/white colony screening.

TOP10 competent cell genotypes: F-*mcrA* Δ (*mrr-hsdRMS-mcrBC*) ϕ 80*lacZ* Δ M15 Δ *lacX74 recA1 ara* Δ 139 Δ (*ara-leu*) 7697 *gaU gaK rpsL* (StrR) *endA1 nupG*.

Components

Components No.	Name	11801ES80	11801ES92
11801	TOP10 Chemically Competent Cell	10 \times 100 μ L	100 \times 100 μ L

Specifications

Species	<i>Escherichia coli</i>
Cell type	Chemical competent cell
Whether it contains an F' episome	No
Efficiency	$>1 \times 10^8$
Blue/white colony screening	Yes
Bacterial or yeast strains	TOP10
Amplification of the ccdB-containing vector	It is not suitable for amplification of ccdB-containing vectors

Shipping and Storage

The product is shipped with dry ice and can be stored at -85°C~-65°C for six months. Do not store the product in -20°C or liquid nitrogen.

Instructions

1. Take 100 μ L of competent cells, ice bath, thaw (approximately 5 min).
2. Add the DNA of interest to the thawed competent cell suspension immediately, gently flick well, and let stand in an ice bath for 25 min.
【Note】 : Do not add more than one-tenth of the volume of competent cell suspension.
3. Place the EP tube in a 42°C water bath for 45 sec, then quickly transfer to ice and let stand for 2 min.
【Note】 : Do not shake during this process, otherwise it will reduce the conversion efficiency.
4. Add about 700 μ L of antibiotic-free LB or 2YT medium to the centrifuge tube, mix well and recover at 37 °C at 200 rpm for 60 min.
5. Collect the bacteria by centrifugation at 5000 rpm for 1 min, leave about 100 μ L of supernatant coating onto a plate containing the corresponding antibiotic, and incubate at 37°C overnight. For blue/white colony screening,



incubate the plate at 37°C for at least 13 h.

Notes

1. Please wear the necessary PPE, such lab coat and gloves, to ensure your health and safety!
2. For research use only.