



ssDNA Assay Kit

Product description

ssDNA Assay Kit is a simple, sensitive and accurate single strand DNA (ssDNA) fluorescence quantitative detection kit with good linear relationship between 1-200 ng. This kit contains fluorescence detection reagent, buffer solution and related ssDNA standards. Before use, dilute the buffer solution of fluorescence detection reagent into working solution, and then add the ssDNA sample to be tested, then use the fluorescence microplate or Qubit Read with a fluorometer. The selectivity of this kit to single stranded DNA is not higher than that to double stranded DNA, but it has good tolerance to conventional pollutants such as proteins, salts, detergents, etc.

Components

Name	Concentration	12645ES60 (100T)	12645ES60 (500T)
ssDNA Reagent	200× in DMSO	250 μL	1250 μL
ssDNA Buffer	Not applicable	50 mL	250 mL
ssDNA Standard 1	0 ng/μL in TE buffer	1 mL	5×1 mL
ssDNA Standard 2	20 ng/μL in TE buffer	1 mL	5×1 mL

Specifications


Assay	ssDNA Quantitation
Excitation/Emission	500/525
For Use With (Equipment)	Qubit Fluorometer, fluorescence microplate
No. of Reactions	100T/500T
Product Line	Qubit Quantitation
Quantitation Range	1-200 ng, 50pg/μL - 200ng/μL
Sample Volume	1 μL to 20 μL
Detection Method	Fluorescent

Shipping and Storage

All the components are shipped with ice packs and can be stored at 2°C~8°C away from light for 6 months. Please avoid repeated freeze-thawing.

Instructions

1. The operation flowchat



2. Preparation

2.1 Equilibrate the components to room temperature before use.

2.2 Prepare a sufficient of 0.5-mL tubes and label the table at the lid. Do not label the side wall of the tube to avoid affecting the fluorescence signal collection. The recommended 0.5-mL tubes include Qubit assay tubes (Cat#Q32856) and Axygen PCR-05-C tubes (Cat#10011-830).

3. Prepare the working solution

Use ssDNA Buffer to dilute the ssDNA Reagent to 1× (for example: take 1 μL of ssDNA Reagent and add 199 μL of ssDNA Buffer). Working solution is prepared for immediately use within 3 h.

4. Prepare the samples to be tested

4.1 Add 190 μL of 1× working solution into 0.5-mL tube, then add 10 μL of ssDNA Standard 1 and ssDNA Standard 2 to the corresponding standard PCR tube. Vortex gently for 2-3 sec and be careful not to create bubbles.

4.2 Add 180-199 μL of 1× working solution and 1-20 μL of each DNA samples to the appropriate 0.5-mL tube. The final volume of each sample in the tube is 200 μL. Vortex gently for 2-3 sec and be careful not to create bubbles.

5. Fluorescence Detection

5.1 Incubate the tube at room temperature in dark for 2 min.

5.2 According to the operating instructions of the Qubit Fluorometer, select ssDNA High Sensitivity to measure the fluorescence signal value.

Notes

For your safety and health, please wear lab coats and disposable gloves for operation.