



## UCF.ME™ UltraNuclease ELISA Kit

### Product description

UltraNuclease, also known as non-restrictive endonuclease, broad-spectrum nuclease, is a non-specific endonuclease derived from *Serratia Marcescens*, which hydrolyzes internal phosphodiester bonds between any nucleotides in nucleic acids to produce 5'-monophosphate oligonucleotides of 2-5 bases in length. It can degrade DNA and RNA in various forms (double-stranded, single-stranded, linear, circular, native or denatured) under a very broad range of conditions (6 M Urea, 0.1 M Guanidine HCl, 0.4% Triton X-100, 0.1% SDS, 1 mM EDTA, 1 mM PMSF) and is widely used to remove nucleic acids from biological products. UltraNuclease can also be removed by corresponding methods subsequently.

This kit uses the principle of double-antibody sandwich enzyme-linked immunosorbent assay (sandwich ELISA) to detect the residues of denatured and non-denatured UltraNuclease. First, coat the microplate with anti-UltraNuclease rabbit polyclonal antibody to form a solid-phase antibody. Second, add UltraNuclease standard and test sample to the solid-phase antibody microplate, then add biotin labeled Anti-UltraNuclease polyclonal antibody, and finally add horseradish peroxidase-labeled streptavidin (SA-HRP) to form an antibody + antigen + antibody-Biotin + SA-HRP complex. Subsequently, TMB substrate was added to the complex to observe color reaction after washing the complex. TMB is converted into blue under the catalysis of HRP enzyme and finally converted into yellow in the presence of acid, and the shade of color is positively correlated with the amount of UltraNuclease in the sample. The detection quantification range of this kit is 0.047-3 ng/mL; the lower detection limit is 23.5 pg/mL.

### Components

Components No.	Name	36701ES59 (96 T)
36701-A	ELISA Microplates	1 Plate
36701-B	Detection Antibody: Biotin-conjugated Rabbit Anti-UltraNuclease Antibodies	35 µL (0.2 mg/mL)
36701-C	Standard: UltraNuclease	1 vial (500 ng/mL)
36701-D	HRP-conjugated Streptavidin	10 µL
36701-E	Dilution Buffer 1	45 mL
36701-F	20× Wash Buffer	50 mL
36701-G	Dilution Buffer 2	30 mL
36701-H	TMB	15 mL
36701-I	Stop Solution	10 mL
36701-J	Sealing Film	5 pieces

### Specification

Sensitivity	23 pg/mL (range 0.047~3 ng/mL)
Assay Time	<4 hours
Assay Principle	Two-site immunoenzymetric assay



Signal Amplification	Biotin-Streptavidin system
Detection Wavelength	450nm

## Shipping and Storage

The UltraNuclease ELISA kit is shipped with ice pack and can be stored at 2°C ~8°C. Unopened product is valid for one year. Once the reagent is opened, it is valid for half a year. Never freeze this product.

Storage of the prepared reagents:

1. The diluted washing buffer and dilution buffer can be stored at 2°C ~8°C for 1 week
2. Standard 36701-C is liquid and stored at 2°C ~8°C for 1 year
3. Prepared detection antibody and stop buffer can be stored at 2°C ~8°C for 1 month

## Instructions

### 1. Reagent preparation

All reagent components and samples to be tested need to be returned to room temperature before use.

1.1 Preparation of 1× Wash Buffer: Equilibrate the concentrated solution to room temperature and dissolve it fully without crystals. Mix well and transfer 50 mL of 20× Wash Buffer into container, and then bring the volume to 1 L with ultrapure water. the specific volume can be prepared according to the amount used each time.

1.2 Preparation of Detection Antibody: Centrifuge at 10,000 r for 20 s before use, and then dilute the detection antibody to a working solution concentration of 0.5 µg/mL with Dilution Buffer 2.

1.3 Preparation of HRP-conjugated Streptavidin: Centrifuge at 10,000 r for 20 s before use, and then dilute the HRP-conjugated Streptavidin to the working concentration by 1:5000 dilution with Dilution buffer 2.

1.4 Preparation of the standard curve: Prepare 8 sterile 1.5 mL centrifuge tubes and label them sequentially according to the standard concentration. Pipette 994 µL of Dilution Buffer 1 to the centrifuge tube marked 3 ng/mL, and pipette 500 µL to the other tubes. Transfer 6 µL of 500 ng/mL UltraNuclease Standard to the centrifuge tube marked 3 ng/mL and mix well, then transfer 500 µL to the centrifuge tube with the next labeled concentration, mix well, and perform a series of 2-fold gradient dilutions of the standard. It can be carried out according to the following table and the initial maximum concentration is 3 ng/mL, the minimum is 0.047 ng/mL. A corresponding standard curve needs to be prepared for each test, and standard curves of different kits and different times cannot be mixed. For sample testing, the required standard volume for each well is 100 µL. Note that the preparation volume should be higher than the required volume to avoid insufficient volume.

Table 1. Preparation of UltraNuclease standard system (detected by microplate reader 0.047-3 ng/mL)\*

Vial	Dilution Buffer 1 (µL)	500 ng/mL UltraNucleaseStandard volume (µL)	Standard final concentration (ng/mL)
A	994	6	3
B	500	500 A	1.5
C	500	500 B	0.75
D	500	500 C	0.375
E	500	500 D	0.188
F	500	500 E	0.094



G	500	500 F	0.047
H	300	0	0=Blank




Figure 1. Simplified flow chart of standard configuration

## 2. Assay procedure

All reagent components and samples to be tested need to be returned to room temperature before use. It is strongly recommended that all standards and test samples be assayed in duplicate.

2.1 Reagent preparation: Prepare reagents, diluted standards and samples.

2.2 Strip Determination: Calculate the strips required for test samples and standards, take strips from the foil bag, put the remaining strips back into the foil bag, seal the bag, and store at low temperature.

2.3 ELISA plate washing : Wash the plate three times with 1× Wash Buffer (300 μL/well), and pat the ELISA plate dry. Washing the plate has a major impact on the test results, make sure that no wash solution remains on the plate.

2.4 Samples incubation: Add standard and test samples, 100 μL/well, ensure that the sample adding is completed within 15 minutes, and incubate at 37°C for 1 h.

2.5 ELISA plate washing: Discard the liquid in the well, wash the plate with 1× Wash Buffer (300 μL/well) five times, and pat the ELISA plate dry.

2.6 Biotin-conjugated detection antibody incubation: Add the detection antibody pre-prepared to the working concentration into the ELISA plate, 100 μL/well, and incubate at 37°C for 1 h.

2.7 ELISA plate washing: Discard the liquid in the wells, wash the plate with 1× Wash Buffer (300 μL/well) five times, and pat the ELISA plate dry.

2.8 HRP-conjugated Streptavidin incubation: Add the conjugated HRP avidin pre-prepared to the working concentration into the ELISA plate, 100 μL/well, and incubate at 37°C for 40 min.

2.9 ELISA plate washing: Discard the liquid in the well, wash the plate 1× Wash Buffer (300 μL/well) five times, and pat the ELISA plate dry.

2.10 Color development: Return the substrate solution TMB to room temperature 10 minutes before use, add the TMB to the ELISA plate, 100 μL/well, and incubate at 37°C for 15 minutes in the dark.

2.11 Stopping: Add stop solution to the plate, 50 μL/well, and shake the plate gently until the color is uniform.

2.12 Reading: Read the absorbance at 450 nm within 20 min.

## 3. Analysis of results

3.1 If the OD value of the sample exceeds the maximum peak of the standard curve, the sample should be diluted and re-measured.

3.2 Standard curve drawing: ELISACalc.exe Regression analysis/ Curve fitting software was recommended to draw the standard Curve (X-standard concentration ng/mL; Y-final OD450 nm) by regression of a conic fitting according to the absorbance of the UltraNuclease standard (subtracting the OD value of the standard blank well namely the final



reading, or not subtracting the OD value of the blank well),. Calculate the nuclease concentration in the sample based on the standard curve and the dilution folds of the sample.

Concentration (ng/mL)	OD450 Value
3	2.359
1.5	1.559
0.750	0.960
0.375	0.585
0.188	0.397
0.093	0.326
0.047	0.242
0.023	0.188
0	0.063




Figure 2. Representative standard curve

【Notes】 :This graph is used for reference only. A standard curve should be generated for each assay.

## Notes

1. The kit needs to be used up within its shelf life. Mixing of related reagents from different batches is prohibited.
2. The kit is designed for detecting the target antigens and samples marked in the instructions only. Other applications need to be designed and verified by the user, and the reliability and accuracy should be evaluated based on the results.
3. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.
4. This product is for research use ONLY!