

Hygromycin B

Product Information

Product Name	Cat#	Size
Hygromycin B	60225ES03	1 g
	60225ES10	10 g

Product Description


Hygromycin B, an aminoglycoside antibiotic synthesized by *Streptomyces hygrosopicus*, inhibits protein synthesis by interfering with 70S ribosomal translocation and inducing misreading of mRNA template. Thus killing prokaryotes (such as bacteria), eukaryotes (such as yeasts, fungi) and higher mammalian eukaryotes.

Hygromycin resistance genes (*hyg* or *hph*) derived from *Escherichia coli* encode hygromycin B phosphotransferase, which detoxifies hygromycin B into a non-biologically active phosphorylated product, making it a very useful selective marker for screening and culturing prokaryotic or eukaryotic cells successfully transfected with hygromycin resistance genes. In addition, due to different modes of action, Hygromycin B is often used in combination with G418 or Blasticidin S for selection of dual-resistant positive cell lines. Hygromycin B can also be used as an antiviral agent because it selectively penetrates into cells with increased membrane permeability due to viral infection and has the effect of inhibiting translation. It can also act as an insect-repellant via mixed with animal feed.

Product Properties

CAS No.	31282-04-9
Molecular formula	C ₂₀ H ₃₇ N ₃ O ₁₃
Molecular weight	527.52 g/mol
Appearance	White to light yellow to brown powder
Purity	>90% (HPLC)
Solubility	soluble in H ₂ O
Potency	≥1050 U/mg

Structure



Shipping and Storage

The product is shipped with ice pack and can be stored at -20°C for 2 years.

Cautions

1. Hygromycin B resistance gene (*hyg* or *hph*), except from *E. Coli.*, are also found in other bacterial strains, including *Streptomyces hygrosopicus* and *Klebsiella pneumoniae*.
2. Toxic compound, avoid skin and eyes contact, please handle with care.
3. For your safety and health, please wear lab coats and disposable gloves for operation.
4. For research use only!

Instructions

设置格式[10784]: 两端对齐, 段落间距段前: 0.5 行, 段后: 0.5 行

1. Preparation of storage solution (50 mg/mL)

Weigh 0.5 g hygromycin B and add 10 mL 1×PBS, PH 7.4. After fully dissolved, filtrate with a 0.22 μm filter for sterilization.

Aliquot the sterilized solution and store at -20°C.

2. Commonly used screening concentration

The working concentration of hygromycin B for screening stable strains varies according to cell type, culture medium, growth conditions and cell metabolic rate. It is recommended to establish a kill curve (dose-response curve), to determine the optimal screening concentration for the first time.

Generally, mammalian cells: 50-500 μg/mL; Bacterial/plant cells: 20-200 μg/mL; Fungi: 300-1000 μg/mL.

3. Establishment of killing curve

【Note】 In order to screen stable cell lines, it is necessary to determine the minimum concentration of antibiotics capable of killing untransfected host cells. This can be achieved by establishing a killing curve (dose-response curve). At least five concentrations should be arranged.

1) Day 1: Untransformed cells are plated in an appropriate culture plate at a cell density of 20-25% and cultured overnight.

【Note】 The amount of inoculation cells can be increased for cells requiring higher density to detect vitality.

2) Set the concentration gradient within the appropriate range according to the cell type. Mammalian cell can set 50, 100, 250, 500, 750, 1000 μg/mL. Dilute the Hygromycin B solution 1:10 with deionized water or PBS buffer to 5 mg/mL. And then dilute the solution to the corresponding working concentration according to the following table.

Final Concentration (μg/mL)	Medium Volume (mL)	Addition Volume of 5 mg/mL Hygromycin B (mL)
50	9.9	0.1
100	9.8	0.2
250	9.5	0.5
500	9.0	1.0
750	8.5	1.5
1000	8.0	2.0

3) Day 2: Replace with a freshly prepared medium containing the corresponding concentration of the drug. Make three parallel samples for each concentration.

4) Replace with fresh media containing drugs every 3-4 days.

5) Living cell counts are performed at a fixed cycle (e.g., every 2 days) to determine the appropriate concentration to prevent the growth of untransfected cells. Choose the minimum concentration which kills the majority of cells within an ideal number of days (usually 7-10 days), as the working concentration for screening stable cell line.

4. Screening of stable transfected cells

1) 48 h after transfection, the cells were subcultured by screening medium containing hygromycin B at appropriate concentration (direct or diluted).

【Note】 Antibiotics work best on actively dividing cells. If the cells are too dense, the antibiotic will not kill the cells. Split the cells such that the cells are no more than 25% confluent.

2) Change the screening medium every 3-4 days.

3) Measure the cell colony-formation after 7 days of screening. Colony formation may take another week or more, depending on host cell type, transfection, and screening effectiveness.

4) Pick 5-10 resistant clones and transfer to 35 mm cell culture plates, and cultured in drug-containing screening medium for 7 days.

5) Replace with fresh medium without drugs for culture.