

Hieff[®] Taq DNA Polymerase

Product Information

Product Name	Cat#	Size
Hieff [®] Taq DNA Polymerase	10101ES80	1,000 U
	10101ES92	10,000 U

Product Description

Hieff[®] Taq DNA polymerase is a thermostable recombinant DNA polymerase expressed by *Thermus aquaticus* with a molecular weight of 94 KDa. It has 5'→3' polymerase activity and 5'→3' exonuclease activity but without 3'→5' exonuclease activity. The amplified product has 3'-dA and can be directly used for TA cloning.

Product Components

Component Number	Name	Cat#/Size	
		10101ES80(1,000 U)	10101ES92(10,000 U)
10101-A	10×Taq Buffer (Mg ²⁺ Free)	4×1 mL	4×10 mL
10101-B	25 mM MgCl ₂	2×1 mL	2×10 mL
10101-C	Hieff [®] Taq DNA Polymerase (5 U/μL)	200 μL	2×1 mL

Application

Genotyping, colony PCR and other conventional PCR.

Unit Definition

Using the activated DNA of salmon sperm as template/primer, the activity is defined as one active unit (U) when 10 nmol of total nucleotide was ingested as acid insoluble substance at 74°C for 30 min.

Shipping and Storage

The products are shipped with ice pack and can be stored at 20°C for 2 years.

Product Notes

- 1 For your safety and health, please wear lab coats and disposable gloves for operation.
- 2 For research use only!

Reaction composition(Preparation on ice)

Components	Volume(μL)	Final Concentration
ddH ₂ O	to 50	-
10×Taq Buffer (Mg ²⁺ Free)	5	1×
25 mM MgCl ₂	3	1.5 mM
dNTP Mix (10 mM each)	1 μL	0.2 mM
DNA template	X μL	-
Forward primer (10 μM)	2 μL	0.4 μM
Reverse primer (10 μM)	2 μL	0.4 μM
Hieff [®] Taq DNA Polymerase (5 U/ μL)	0.4 μL	0.04 U/ μL

Notes:

- 1) Final concentration of Mg²⁺: The optimal concentration of Mg²⁺ is 1.5-2 mM. If necessary, the optimal concentration of Mg²⁺ can be explored upward at intervals of 0.2-0.5 mM.
- 2) Polymerase addition: The polymerase has a certain degree of 5'- 3' polymerase activity at room temperature. In order to prevent non-specific amplification, it is suggested to add the polymerase to the reaction system in the last step.
- 3) Concentration of polymerase: 0.04 U/ μL is recommended. It can be optimized between 0.025-0.04 U/ μL .
- 4) Recommended use of different templates (50 μL reaction system)

Type of template	Template usage
Genomic DNA	50 ng-100 ng
Plasmid DNA	10 pg-20 ng
cDNA	1-5 μL (No more than 1/10 of the reaction system)

Thermal cycling protocol

Stage	Temperature (°C)	Time	Cycles
Pre-denaturation	94	30 sec-5 min	1
Denaturation	94	30 sec	} 35
Annealing	50-60	30 sec	
Extension	72	60 sec/kb	
Final Extension	72	10 min	1

Notes:

- 1) Pre-denaturation temperature and time: 94°C is recommended. The recommended pre denaturation time: 30 sec for plasmid DNA and other simple templates; 3 min for complex templates such as cDNA and genomic DNA; 5-10 min for the template with high GC.
- 2) Annealing temperature and time: 60°C is recommended. Temperature gradient can be set up to find the optimum temperature for primer annealing. The recommended annealing time is set to 20 sec and can be adjusted within 10-30 sec. Too long annealing time may cause the amplified products diffusion on the agarose gel.
- 3) Amplification products: Please store the PCR amplification products at - 20°C to prevent DNA degradation.