


UCF.ME® UltraNuclease, GMP-grade

About YEASEN



What is UCF.ME® UltraNuclease?

UltraNuclease is a nonspecific endonuclease that cleavages and digests DNA and RNA of all formats (double-stranded, single-stranded, linear, circular, naturally occurring, or denatured) under a broad range of experimental conditions. It has been widely used for nucleic acid removal for numerous applications.







Vaccine Manufacturing



Drug Purification



Cell Clumping Prevention

With our dedicated R&D team's over 3 years of continuous refining, Yeasen has developed the gene-modified E. Coli-expressed UCF.ME® UltraNuclease, which is characterized by its exceptional purity and high nuclease activity. Manufactured under the GMP guideline, this product has no protease activity and is free of endotoxin or viral contaminants. It can be applied not only for scientific research but also for nucleic acid removal during biological products' manufacturing process.


6 Reasons for Choosing UCF. ME® UltraNuclease

- ▶ 1.Digests DNA and RNA of all formats
- ▶ 2.Ultra-high nuclease activity; Reduces nucleic acid residues to pg-level within short periods of treatment.
- > 3. Highly thermostable and stress-tolerant, compatible with various experimental conditions.
- ▶ 4.Non-animal origin, free of probiotics or endotoxins.
- ▶ 5.Manufactured under GMP guidelines; Capable of both research and scaled-up manufacturing.
- ▶ 6.Complete audit documentation, enabling smooth product certification application.



UCF.ME® UltraNuclease Nuclease Activity Quantification

▶ UltraNuclease's specific activity quantification using universal substrate (>1.5×106U/mg)



* The definition of one Activity unit (U) is the amount of enzyme used to change the absorption value of \triangle A260 by 1.0 in 30 minutes in a 2.625 mL reaction system at 37°C with a pH of 8.0 (equivalent to complete digestion of 37 µg salmon sperm DNA into oligonucleotides).


UCF.ME® UltraNuclease Purity Test

▶ SDS-PAGE purity test result: ≥ 99%

kDa 180 140 100 75 60 45 35 25 15

SDS-PAGE purity test result


▶ SEC-HPLC purity test result: ≥ 99%



Peak	Retention Time (min)	Peak Area (%)
Minor Peak	18.328	0.79
Major Peak	19.850	99.21

SEC-HPLC purity test result

Application and Post-treatment Residues



Cell Lysate Viscosity Reduction
Tube A: Treated with 25U/mL UltraNuclease
Tube B: Without UltraNuclease Treatment

UCF.ME® UltraNuclease Quality Inspection

Catalog Number	20156ES	20157ES
Manufacturing Grade	RUO-grade	GMP-grade
Expression Host	E. coli	E. coli
Animal Originicity	yes	no
Antibiotics	yes	no
IPTG Residue Verification	_	\checkmark
Manufacturing Standard	_	GMP
Appearance (Colorless)	√	\checkmark
Specific Activity	≥1.5×10 ⁶ U/mg	≥1.5×10 ⁶ U/mg
Protein Content (BCA)	≥0.1mg/mL	≥0.1mg/mL
Protein Purity SDS-PAGE (≥99%)	√	\checkmark
Protein Purity SEC-UPLC (≥99%)	_	\checkmark
Proteinase Residues (universal substrate)	_	✓
Bacterial Endotoxin (gel-clot)	_	<0.25EU/1000U
Host Protein Residues (ELISA)	_	<10ppm (µg/mL)
1101 sterility test	_	\checkmark
Pathogen (HIV/HBV/HCV) Detection (PCR)	_	✓

Tested and Qualified (\checkmark); Test not performed (-)

Product Ordering Information

Product Name	Catalog Number	Size
UCF.ME® UltraNuclease	20156ES25/50/60	25/50/100KU
UCF.ME® UltraNuclease, GMP-grade	20157ES25/50/60	25/50/100KU

Visit our official website: www.yeasen.com for more information...

Yeasen Biotechnology (Shanghai) Co., Ltd.

Room 202, Building No. 2, Lane 166, Tianxiong Road, Shanghai International Medical Zone, Pudong New Area, Shanghai, China

Service Hotline: 400-6111-883 Web: www.yeasen.com Technical Support: techserv@yeasen.com
Product consultation: marketing@yeasen.com



