

10× Transcription Buffer GMP-grade (Mg²⁺ free)

Product Information

Product Name	Cat#	Size
	10669ES03	1 mL
10×Transcription Buffer GMP-grade (Mg ²⁺ free)	10669ES10	10 mL
	10669ES60	100 mL

Product Description

This product is optimized transcription buffer and efficient for in vitro synthesis of RNA. The formulation of the product: 400mM Tris-HCl, 20mM spermidine, 100mM DTT, pH 7.9 at 25°C.

This product is produced in accordance with GMP regulations and provided in liquid form.

Contents

Contents No.	Name	Catalog No./Specification		
		10669ES03	10669ES10	10669ES60
10669	10×Transcription Buffer GMP-grade (Mg ²⁺ free)	1 mL	10 mL	100 mL

Shipping and Storage

This product is shipped with dry ice and can be stored at -15°C ~ -25°C for one year.

Notes

1. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.