

10× Transcription Buffer GMP-grade

Product Information

Product Name	Cat#	Size
10×Transcription Buffer GMP-grade	10627ES03	1 mL
	10627ES10	10 mL
	10627ES25	25 mL
	10627ES76	500 mL

Product Description

This product is optimized transcription buffer and efficient for in vitro synthesis of RNA. The formulation of the product: 400 mM Tris-HCl, 460 mM MgCl₂, 20 mM spermidine, 100 mM DTT, pH 7.9 at 25°C.

This product is produced in accordance with GMP regulations and provided in liquid form.

Contents

Contents No.	Name	Catalog No./Specification			
		10627ES03	10627ES10	10627ES25	10627ES76
10627	10× Transcription Buffer GMP-grade	1 mL	10 mL	25 mL	500 mL

Shipping and Storage

This product is shipped with dry ice and can be stored at -15°C ~ -25°C for one year.

Notes

1. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.