



Ver. HB221111

10× Capping Buffer GMP-grade

Product description

This product is optimized 10× capping buffer and efficient for mRNA capping reaction and 5'end labeling reaction.

The formulation of the product: 0.5 M Tris-HCl, 50 mM KCl, 10 mM MgCl₂, 10 mM DTT, pH 8.0 at 25°C

This product is produced in accordance with GMP regulations and provided in a liquid form.

Components

Components No.	Name	10666ES03	10666ES10	10666ES25	10666ES76
10666	10× Capping Buffer GMP-grade	1 mL	10 mL	25 mL	500 mL

Shipping and Storage

This product is shipped with dry ice and can be stored at -15°C ~ -25°C for one year.

Notes

1. For your safety and health, please wear personal protective equipment (PPE), such as laboratory coats and disposable gloves, when operating with this product.