

HB211208

L-alanyl-L-glutamine solution, 200mM

Product Information

Product name	Cat#	Size
L-alanyl-L-glutamine solution,200mM	60701ES20	20 mL
L-alanyl-L-glutamine solution,200mM	60701ES60	100 mL
L-alanyl-L-glutamine solution,200mM	60701ES76	500 mL

Product Description

L-alanyl-L-glutamine is a substitute for L-glutamine, which is widely used in cell culture. Due to the high stability in the aqueous solution, it can maintain a long-term and stable supply of nutrients and reduce the cytotoxicity caused by decomposition products. L-alanyl-L-glutamine contains 0.85% (W:V) NaCl, which can be used as an additive for mammalian adherent cells or suspension cell culture.

Product Properties

Product name	L-alanyl-L-glutamine solution, 200mM
Appearance	liquid
Molecular Formula	$C_8H_{15}N_3O_4$
Molecular Weight	217.22

Shipping and Storage

The product is shipped with ice pack and can be stored at 2-8°C for two years.

Caution

For your safety and health, please wear lab coats and disposable gloves for operation. For research use only!